

Product data sheet

Rev 00 / 27.04.2022

315130-000

Product name EN

HP 120 RA patient line, PVC-flex

Product group

4 - Patientenschläuche 83 bar



FlowRate
0,1-30 ml/s

Specifications

Highlights

**83 bar/
1200 psi**
max. Pressure

**Accutron HP
Accutron HP-D**
Compatible Devices

3 years
Shelf life

Intended use:

The patient lines are intended for the fluid transfer of contrast medium (CM) and saline solution (NaCl) into the patient in angiography.

Application and Application frequency:

Single use device

Single use device intended for single use on a single patient.

Net weight (Product + Blister): 36,3 g

Carton weight: 295 g

Quantity per outer carton: PU 50

Minimum order quantity = 1 full carton/ box: PU 50


Weight Shipping unit: // kg

Product data sheet

Rev 00 / 27.04.2022

Technical Specification

Functional facts



Graphic

Technical Specification

Material data

Material Tube
 Tube: PVC, no DEHP
 Luer-Lock connector: ABS, MABS
 Screw cap: PP
 Protective cap: white
 Rotating adapter: PC

Material Syringe
 NA

Material Accessories
 NA

Sterile barrier system
 Peel off soft blister packaging

Materials sterile barrier system
 Special medical paper 60g/m²
 PA/PE multilayer coextruded film

Dimensions outer carton
 38,9x34x11,2 cm

CE-mark
 CE 0482

Classification (MDD)
 Class IIa

UMDNS Code
 15-265

GMDN Code
 16545

Sterilization
 ETO-Gas

Storage conditions
 Storage temperature: 10°C – 30°C
 Relative humidity: 30% r.L. – 80% r.L.
 Store in the closed carton,
 protected from direct sunlight and
 in a dry place

HS Code
 90183900

Barcode inner blister
 (01)14250288600307(10)1234567
 89(17)yymmdd

Barcode inner carton
 NA

Barcode outer carton
 (01)34250288600301(10)1234567
 89(17)yymmdd

Manufacturer
MEDTRON AG
Hauptstrasse 255
66128 SAARBRÜCKEN
Germany